

# Tung BRUSH™

**Eliminates bad breath  
where it starts.**

Clinical research shows that **odor-causing bacteria** on the tongue is the **#1** cause of bad breath. There's no better tool to get rid of that nasty bacteria and "**tongue gunk**" than the Tung Brush™.



**LICKS BAD BREATH  
& TONGUE GUNK.**

Designed to outperform tooth brushes and scrapers with **greater coverage** and **superior comfort** with no gagging.

Unique bristle design penetrates tongue crevices not visible to the naked eye where odor-causing bacteria hide.

[www.tungbrush.com](http://www.tungbrush.com)

**Directions:** Extend tongue and gently brush from the back of the tongue forward to the tip (for approx. 10 seconds). Repeat as necessary. Replace your Tung Brush™ at least every three months.

Should not be used on teeth and gums. Slight tenderness and/or bleeding may occur on tongue at first, but should stop with regular brushing. If either condition continues, stop use and consult your dentist. Excessive brushing may cause soreness.

*TungGel™, TungBrush™ and the logos are trademarks of Peak Enterprises, Inc.  
All rights reserved.*

Made in the U.S.A.  
U.S. Patents RE37,625 / 404,129  
Copyright 2006 Peak Enterprises, Inc.  
Peak Enterprises, Inc. Sarasota, FL USA

